

## Product datasheet for **RC210190L4V**

### UBASH3A (NM\_018961) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | UBASH3A (NM_018961) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | UBASH3A  |
| Synonyms:                 | CLIP4; STS-2; TULA; TULA-1   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_018961  |
| ORF Size:                 | 1983 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC210190).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_018961.2</a>  |
| RefSeq Size:              | 2492 bp  |
| RefSeq ORF:               | 1986 bp  |
| Locus ID:                 | 53347  |
| UniProt ID:               | <a href="#">P57075</a>   |
| Cytogenetics:             | 21q22.3  |
| Domains:                  | UBA, PGAM, SH3   |
| Protein Families:         | Druggable Genome   |



[View online »](#)

**MW:** 73.9 kDa

**Gene Summary:** This gene encodes one of two family members belonging to the T-cell ubiquitin ligand (TULA) family. Both family members can negatively regulate T-cell signaling. This family member can facilitate growth factor withdrawal-induced apoptosis in T cells, which may occur via its interaction with AIF, an apoptosis-inducing factor. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Aug 2011]