

Product datasheet for **RC210214L4V**

SLC7A2 (NM_001008539) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC7A2 (NM_001008539) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC7A2
Synonyms:	ATRC2; CAT2; HCAT2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008539
ORF Size:	1974 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210214).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008539.1
RefSeq Size:	7389 bp
RefSeq ORF:	1977 bp
Locus ID:	6542
UniProt ID:	P52569
Cytogenetics:	8p22
Protein Families:	Druggable Genome, Transmembrane
MW:	71.5 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a cationic amino acid transporter and a member of the APC (amino acid-polyamine-organocation) family of transporters. The encoded membrane protein is responsible for the cellular uptake of arginine, lysine and ornithine. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]