

Product datasheet for **RC210224L3V**

Cannabinoid Receptor II (CNR2) (NM_001841) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cannabinoid Receptor II (CNR2) (NM_001841) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cannabinoid Receptor II
Synonyms:	CB-2; CB2; CX5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001841
ORF Size:	1080 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210224).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001841.1
RefSeq Size:	1776 bp
RefSeq ORF:	1083 bp
Locus ID:	1269
UniProt ID:	P34972
Cytogenetics:	1p36.11
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction




[View online »](#)

MW: 39.5 kDa

Gene Summary: The cannabinoid delta-9-tetrahydrocannabinol is the principal psychoactive ingredient of marijuana. The proteins encoded by this gene and the cannabinoid receptor 1 (brain) (CNR1) gene have the characteristics of a guanine nucleotide-binding protein (G-protein)-coupled receptor for cannabinoids. They inhibit adenylate cyclase activity in a dose-dependent, stereoselective, and pertussis toxin-sensitive manner. These proteins have been found to be involved in the cannabinoid-induced CNS effects (including alterations in mood and cognition) experienced by users of marijuana. The cannabinoid receptors are members of family 1 of the G-protein-coupled receptors. [provided by RefSeq, Jul 2008]

Product images:



[RC210224L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC210224L3V particle to overexpress human CNR2-Myc-DDK fusion protein.