

Product datasheet for **RC210238A1V**

Human Fetuin B (FETUB) (NM_014375) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Fetuin B (FETUB) (NM_014375) AAV Particle
Tag:	Myc-DDK
Symbol:	Fetuin B
Synonyms:	16G2; Gugu; IRL685
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC210238 representing NM_014375 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGGTCTGCTCCTTCCCCTGGCACTCTGCATCCTAGTCCTGTGCTGCGGAGCAATGTCTCCACCCAGC
TGGCCCTCAACCCCTCGGCTCTGCTCTCCCGGGCTGCAATGACTCAGATGTGCTGGCAGTTGCAGGCTT
TGCCCTGCGGGATATTAACAAAGACAGAAAGGATGGCTATGTGCTGAGACTCAACCGAGTGAACGACGCC
CAGGAATACAGACGGGGTGGCTGGGATCTCTGTTCTATCTTACACTGGATGTGCTAGAGACTGACTGCC
ATGTGCTCAGAAAGAAGGCATGGCAAGACTGTGGAATGAGGATATTTTTGAATCAGTTTATGGTCAATG
CAAAGCAATATTTTATATGAACAACCCAGTAGAGTTCTCTATTTAGCTGCTTATAACTGTACTCTTCGC
CCAGTTTCAAAAAAAGATTTACATGACGTGCCCTGACTGCCCAAGCTCCATACCCACTGACTCTTCCA
ATCACCAAGTGCTGGAGGCTGCCACCGAGTCTCTTGCAGAAATACAACATGAGAACACATCCAAGCAGTA
TTCTCTTTCAAAGTCACCGAGGCTTCTAGCCAGTGGGTGGTGGCCCTTCTTACTTTGTGGAATACTTA
ATTAAGAATCACCATGTACTAAATCCCAGGCCAGCAGCTGTTCACTTCAGTCTCCGACTCTGTGCTG
TTGGTCTTTGAAAGTTCTCTGACTCGAACACACTGGGAAAAGTTGTCTCTGTGACTTGTGACTTCTT
TGAATCACAGGCTCCAGCCACTGGAAGTGAAGTCTGCTGTTAACCAGAAACCTACAAACCTTCCCAAG
GTGGAAGAATCCAGCAGAAAAACCCCCCAACAGACTCCCCCTCAAAGCTGGGCCAAGAGGATCTG
TCCAATATCTTCTGACTTGGATGATAAAATTCAGGAAAAGGGCCCTCAGGAGGCCCTTCTGTGCA
TCTGGACCTAACACGAATCCCAGGGAGAAACCTGGATATTTCTTCTCTTCTTCTGGAGCCTATGGAG
GAGAAGCTGGTGGTCTGCCTTTCCCAAGAAAAAGCACGCACTGCTGAGTGCCAGGGCCAGCCAGCA
ATGCCAGCCCTTGTCTTCCGCCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210238 representing NM_014375 Red=Cloning site Green=Tags(s) MGLLLPLALCILVLCGGAMSPQLALNPSALLSRGCNDSVDLAVAGFALRDINKDRKDGYVLRNLNRVND QEYRRGGLGSLFYLTLDVLETCHVLRKKAWQDCGMRIFVESVYGQCKAIFYMNNPSRVLYLAAYNCTLR PVSKKKIYMTCPDCPSSIPTDSSNHQVLEAATESLAKYNNENTSKQYSLFKVTRASSQWVVGPSYFVEYL IKESPCTKSQASSCSLQSSDSVPVGLCKGSLTRTHWEKFSVTCDFEFESQAPATGSENSAVNQKPTNLPK VEESQKQNTPTDPSKAGPRGSVQYLPDLDDKNSQEKGPQEAFPVHLDLTTNPQGETLDISFLFLEPME EKLVLVLPFPKEKARTAECPGAQNASPLVLPP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_014375
ORF Size:	1146 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_014375.2
RefSeq Size:	1640 bp
RefSeq ORF:	1149 bp
Locus ID:	26998
UniProt ID:	Q9UGM5
Cytogenetics:	3q27.3
MW:	41.9 kDa