

Product datasheet for **RC210243L2V**

NANOG (NM_024865) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | NANOG (NM_024865) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | NANOG |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_024865 |
| ORF Size: | 915 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC210243). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_024865.1 |
| RefSeq Size: | 2114 bp |
| RefSeq ORF: | 918 bp |
| Locus ID: | 79923 |
| UniProt ID: | Q9H9S0 |
| Cytogenetics: | 12p13.31 |
| Protein Families: | Cancer stem cells, Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Stem cell - Pluripotency |
| MW: | 34.4 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is a DNA binding homeobox transcription factor involved in embryonic stem (ES) cell proliferation, renewal, and pluripotency. The encoded protein can block ES cell differentiation and can also autorepress its own expression in differentiating cells. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]