

## Product datasheet for **RC210266L4V**

### **RNPEPL1 (NM\_018226) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RNPEPL1 (NM_018226) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNPEPL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018226
ORF Size:	1482 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210266).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018226.3</a> , <a href="#">NP_060696.3</a>
RefSeq Size:	3048 bp
RefSeq ORF:	2178 bp
Locus ID:	57140
UniProt ID:	<a href="#">Q9HAU8</a>
Cytogenetics:	2q37.3
Domains:	Peptidase_M1
Protein Families:	Druggable Genome
MW:	55.5 kDa


[View online »](#)

**Gene Summary:**

Broad specificity aminopeptidase which preferentially hydrolyzes an N-terminal methionine, citrulline or glutamine.[UniProtKB/Swiss-Prot Function]