

Product datasheet for **RC210268L3V**

KLC4 (NM_201521) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KLC4 (NM_201521) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLC4
Synonyms:	bA387M24.3; KNSL8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_201521
ORF Size:	1857 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210268).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_201521.1
RefSeq Size:	2404 bp
RefSeq ORF:	1860 bp
Locus ID:	89953
UniProt ID:	Q9NSK0
Cytogenetics:	6p21.1
Protein Families:	Druggable Genome
MW:	68.6 kDa


[View online »](#)

Gene Summary:

Kinesin is a microtubule-associated force-producing protein that may play a role in organelle transport. The light chain may function in coupling of cargo to the heavy chain or in the modulation of its ATPase activity (By similarity).[UniProtKB/Swiss-Prot Function]