

Product datasheet for **RC210274L3V**

Caveolin 1 (CAV1) (NM_001753) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Caveolin 1 (CAV1) (NM_001753) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Caveolin 1
Synonyms:	BSCL3; CGL3; LCCNS; MSTP085; PPH3; VIP21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001753
ORF Size:	534 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210274).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001753.3
RefSeq Size:	2704 bp
RefSeq ORF:	537 bp
Locus ID:	857
UniProt ID:	Q03135
Cytogenetics:	7q31.2
Domains:	Caveolin
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: Focal adhesion, Viral myocarditis

MW: 20.3 kDa

Gene Summary: The scaffolding protein encoded by this gene is the main component of the caveolae plasma membranes found in most cell types. The protein links integrin subunits to the tyrosine kinase FYN, an initiating step in coupling integrins to the Ras-ERK pathway and promoting cell cycle progression. The gene is a tumor suppressor gene candidate and a negative regulator of the Ras-p42/44 mitogen-activated kinase cascade. Caveolin 1 and caveolin 2 are located next to each other on chromosome 7 and express colocalizing proteins that form a stable hetero-oligomeric complex. Mutations in this gene have been associated with Berardinelli-Seip congenital lipodystrophy. Alternatively spliced transcripts encode alpha and beta isoforms of caveolin 1.[provided by RefSeq, Mar 2010]