

## Product datasheet for **RC210279L4V**

### MARK2 (NM\_004954) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MARK2 (NM_004954) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MARK2
Synonyms:	EMK-1; EMK1; PAR-1; Par-1b; Par1b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004954
ORF Size:	2172 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210279).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004954.3</a>
RefSeq Size:	4556 bp
RefSeq ORF:	2175 bp
Locus ID:	2011
UniProt ID:	<a href="#">Q7KZ17</a>
Cytogenetics:	11q13.1
Domains:	UBA, pkinase, KA1
Protein Families:	Druggable Genome, Protein Kinase


[View online »](#)

**MW:** 81.2 kDa

**Gene Summary:** This gene encodes a member of the Par-1 family of serine/threonine protein kinases. The protein is an important regulator of cell polarity in epithelial and neuronal cells, and also controls the stability of microtubules through phosphorylation and inactivation of several microtubule-associating proteins. The protein localizes to cell membranes. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2009]