

Product datasheet for RC210290A1V

Human SNRPB (NM_003091) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNRPB (NM_003091) AAV Particle
Tag: Myc-DDK
Symbol: SNRPB
Synonyms: CCMS; COD; Sm-B/B'; SmB/B'; SmB/SmB'; snRNP-B; SNRPB1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210290 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACGGTGGGCAAGAGCAGCAAGATGCTGCAGCATATTGATTACAGGATGAGGTGCATCCTGCAGGACG
 GCCGGATCTTCATTGGCACCTTCAAGGCTTTTGACAAGCACATGAATTTGATCCTCTGTGACTGTGATGA
 GTTCAGAAAGATCAAGCCAAAGAACTCCAAACAAGCAGAAAGGGAAGAGAAGCGAGTCCTCGGTCTGGTG
 CTGCTGCGAGGGGAGAATCTGGTCTCAATGACAGTAGAGGGACCTCCTCCCAAAGATACTGGTATTGCTC
 GAGTTCACCTTCTGGAGCTGCCGGGGGCCAGGGATCGGCAGGGCTGCTGGCAGAGGAATCCCACTGG
 GGTTCCCATGCCCCAGGCTCCTGCAGGACTTCTGCGCCAGTCCGTGGGGTTGGCGGGCCATCCCAACAG
 GTGATGACCCCAAGGAAGAGTACTGTTGCAGCCGCTGCAGCTGCTGCCACAGCCAGTATTGCCGGGG
 CTCCAACCCAGTACCCACCTGGCCGTGGGGTCTCCCCACCTATGGCCGAGGAGCACCCCTCCAGG
 CATGATGGGCCACCTCCTGGTATGAGACCTCTATGGGTCCCCAATGGGGATCCCCCTGGAAGAGGG
 ACTCCAATGGGCATGCCCTCCGGGAATGCGGCCTCCTCCCCCTGGGATGCGAGGCCTTCTT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210290 protein sequence
 Red=Cloning site Green=Tags(s)

MTVGKSSKMLQHIDYRMRCILQDGRIFIGTFKAFDKHMLILCDCDEFKIKPKNSKQAEREERVLGLV
 LLRGENLVSMTVEGPPPKDTGIARVPLAGAAGGPGIGRAAGRGIPAGVPMPQAPAGLAGPVRGVGGPSQQ
 VMTPQGRGTVAATAAATASIAGAPTQYPPGRGGPPPMGRGAPPPGMMGPPPMGRPPPMGPPMIPPGRG
 TPMGMPPPMRPPPPGMRGLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_003091
ORF Size:	693 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003091.3
RefSeq Size:	1153 bp
RefSeq ORF:	696 bp
Locus ID:	6628
UniProt ID:	P14678
Cytogenetics:	20p13
MW:	23.7 kDa