

Product datasheet for **RC210292L1V**

GNA11 (NM_002067) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GNA11 (NM_002067) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GNA11
Synonyms:	FBH; FBH2; FHH2; GNA-11; HHC2; HYPOC2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_002067
ORF Size:	1077 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210292).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002067.1
RefSeq Size:	4145 bp
RefSeq ORF:	1080 bp
Locus ID:	2767
UniProt ID:	P29992
Cytogenetics:	19p13.3
Protein Pathways:	Calcium signaling pathway, Gap junction, GnRH signaling pathway, Long-term depression, Vascular smooth muscle contraction




[View online »](#)

MW: 42.1 kDa

Gene Summary: The protein encoded by this gene belongs to the family of guanine nucleotide-binding proteins (G proteins), which function as modulators or transducers in various transmembrane signaling systems. G proteins are composed of 3 units: alpha, beta and gamma. This gene encodes one of the alpha subunits (subunit alpha-11). Mutations in this gene have been associated with hypocalciuric hypercalcemia type II (HHC2) and hypocalcemia dominant 2 (HYPOC2). Patients with HHC2 and HYPOC2 exhibit decreased or increased sensitivity, respectively, to changes in extracellular calcium concentrations. [provided by RefSeq, Dec 2013]

Product images:



[RC210292L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC210292L1V particle to overexpress human GNA11-Myc-DDK fusion protein.