

Product datasheet for **RC210294L1V**

SHARP2 (BHLHE40) (NM_003670) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SHARP2 (BHLHE40) (NM_003670) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SHARP2 |
| Synonyms: | BHLHB2; Clast5; DEC1; HLHB2; SHARP-2; SHARP2; STRA13; Stra14 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_003670 |
| ORF Size: | 1236 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC210294). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003670.1 |
| RefSeq Size: | 3061 bp |
| RefSeq ORF: | 1239 bp |
| Locus ID: | 8553 |
| UniProt ID: | O14503 |
| Cytogenetics: | 3p26.1 |
| Domains: | HLH, ORANGE |
| Protein Families: | Transcription Factors |



[View online »](#)

Protein Pathways: Circadian rhythm - mammal

MW: 45.5 kDa

Gene Summary: This gene encodes a basic helix-loop-helix protein expressed in various tissues. The encoded protein can interact with ARNTL or compete for E-box binding sites in the promoter of PER1 and repress CLOCK/ARNTL's transactivation of PER1. This gene is believed to be involved in the control of circadian rhythm and cell differentiation. [provided by RefSeq, Feb 2014]