

## Product datasheet for **RC210301L4V**

### Alpha Taxilin (TXLNA) (NM\_175852) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Alpha Taxilin (TXLNA) (NM_175852) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Alpha Taxilin
Synonyms:	IL14; TXLN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175852
ORF Size:	1638 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210301).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175852.3</a>
RefSeq Size:	4843 bp
RefSeq ORF:	1641 bp
Locus ID:	200081
UniProt ID:	<a href="#">P40222</a>
Cytogenetics:	1p35.2
Protein Families:	Druggable Genome
MW:	61.8 kDa



[View online »](#)

**Gene Summary:**

May be involved in intracellular vesicle traffic and potentially in calcium-dependent exocytosis in neuroendocrine cells.[UniProtKB/Swiss-Prot Function]