

Product datasheet for **RC210310L4V**

IL4I1 (NM_152899) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | IL4I1 (NM_152899) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | IL4I1 |
| Synonyms: | FIG1; hIL4I1; LAAO; LAO |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_152899 |
| ORF Size: | 1701 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC210310). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_152899.1 |
| RefSeq Size: | 1798 bp |
| RefSeq ORF: | 1704 bp |
| Locus ID: | 259307 |
| UniProt ID: | Q96RQ9 |
| Cytogenetics: | 19q13.33 |
| Domains: | Amino_oxidase |
| Protein Families: | Druggable Genome |



[View online »](#)

| | |
|--------------------------|--|
| Protein Pathways: | Alanine, aspartate and glutamate metabolism, Cysteine and methionine metabolism, Metabolic pathways, Phenylalanine, tyrosine and tryptophan biosynthesis, Phenylalanine metabolism, Tryptophan metabolism, Tyrosine metabolism, Valine, leucine and isoleucine degradation |
| MW: | 60.7 kDa |
| Gene Summary: | This gene encodes a secreted L-amino acid oxidase protein which primarily catabolizes L-phenylalanine and, to a lesser extent, L-arginine. The expression of this gene is induced by the cytokine interleukin 4 in B cells. This gene is also expressed in macrophages and dendritic cells. This protein may play a role immune system escape as it is expressed in tumor-associated macrophages and suppresses T-cell responses. This protein also contains domains thought to be involved in the binding of flavin adenine dinucleotide (FAD) cofactor. Multiple transcript variants encoding different isoforms have been found for this gene. Some transcripts of this gene share a promoter and exons of the 5' UTR with the overlapping NUP62 gene. [provided by RefSeq, Jul 2020] |