

Product datasheet for **RC210313L2V**

Ephrin A3 (EFNA3) (NM_004952) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ephrin A3 (EFNA3) (NM_004952) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Ephrin A3
Synonyms:	EFL2; Ehk1-L; EPLG3; LERK3
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004952
ORF Size:	714 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210313).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004952.3
RefSeq Size:	1782 bp
RefSeq ORF:	717 bp
Locus ID:	1944
UniProt ID:	P52797
Cytogenetics:	1q21.3
Domains:	Ephrin
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Axon guidance

MW: 26.4 kDa

Gene Summary: This gene encodes a member of the ephrin (EPH) family. The ephrins and EPH-related receptors comprise the largest subfamily of receptor protein-tyrosine kinases and have been implicated in mediating developmental events, especially in the nervous system and in erythropoiesis. Based on their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are anchored to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. This gene encodes an EFNA class ephrin. [provided by RefSeq, Jul 2008]