

Product datasheet for **RC210321L3V**

OGFRL1 (NM_024576) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	OGFRL1 (NM_024576) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OGFRL1
Synonyms:	dj331H24.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024576
ORF Size:	1353 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210321).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024576.3
RefSeq Size:	1755 bp
RefSeq ORF:	1356 bp
Locus ID:	79627
UniProt ID:	Q5TC84
Cytogenetics:	6q13
Domains:	OGFr_N
Protein Families:	Druggable Genome



[View online »](#)

MW: 51.1 kDa