

## Product datasheet for **RC210335L4V**

### eIF1A (EIF1AX) (NM\_001412) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | eIF1A (EIF1AX) (NM_001412) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | eIF1A  |
| Synonyms:                 | eIF-1A; eIF-4C; EIF1A; EIF1AP1; EIF4C  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001412  |
| ORF Size:                 | 432 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC210335).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001412.3</a>  |
| RefSeq Size:              | 4431 bp  |
| RefSeq ORF:               | 435 bp   |
| Locus ID:                 | 1964   |
| UniProt ID:               | <a href="#">P47813</a>   |
| Cytogenetics:             | Xp22.12  |
| Domains:                  | eIF-1a   |
| MW:                       | 16.3 kDa   |



[View online »](#)

**Gene Summary:**

This gene encodes an essential eukaryotic translation initiation factor. The protein is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNA<sub>i</sub> and eIF3) to the 5' end of capped RNA. [provided by RefSeq, Jul 2008]