

Product datasheet for **RC210358A1V**

Human ART4 (NM_021071) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ART4 (NM_021071) AAV Particle
Tag:	Myc-DDK
Symbol:	ART4
Synonyms:	ARTC4; CD297; DO; DO/ART4; DOK1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC210358 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGTCCATTGATCAACAGATGCAAGAAGATTCTTCTCCAACTACTGTACCTCCTGCAACGATGAGAA
TCTGGCTCCTTGGAGGCTGCTGCCATTCCTGCTGCTCCTCTCTGGCCTGCAGAGACCCACAGAGGGTTC
TGAGGTTGCAATTAATAATCGACTTCGACTTCGCACCAGGTTCTTTGATGATCAGTACCAAGGCTGTAGC
AAACAGGTTATGGAGAACTAACTCAAGGGGATTATTTACAAAAGACATAGAAGCCAGAAGAATTATT
TTAGGATGTGGCAAAAAGCCCACTTAGCCTGGCTTAACCAAGGAAAAGTTCTACCCAGAACATGACTAC
CACACAGCTGTGGCTATTTTGTGTTTACACATTGAACAGCAATGTTCTTCTGACTTTACTAGAGCCATG
GCCTCTGTTGCCAGGACTCCACAGCAGTATGAACGTTTCACTTCAAATATTTACACTACTACCTCA
CCTCAGCAATCCAGCTGCTGAGGAAAGACAGCATCATGGAGAATGGCACTCTGTGCTATGAGGTGCATTA
TAGGACGAAGGATGTCCACTTAAATGCCTACACAGGGGCCACCATTCGATTTGGCCAATTCCTTTCCACA
TCCCTCCTGAAAGAAGAGGCACAGGAGTTTGGGAACCAGACACTATTTACCATATTCACCTGCCTGGGTG
CACCTGTACGTAATCTCCCTCAAGAAGGAAGTCTTGATCCCTCCCTATGAGCTGTTTAAAGTTATAAA
TATGAGCTACCACCAAGAGGAACTGGTTGCAGTTGAGGTCACTGGGAACCTGAGCACATATAACTGT
CAGCTGCTAAAAGCTTCCAGCAAGAAATGCATCCCTGATCCTATAGCTATTGCATCTCTCTCTTTTGA
CCAGTGTCATCATCTTTTCCAAAAGCAGAGTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210358 protein sequence Red=Cloning site Green=Tags(s) MGPLINRCKKILLPTTVPPATMRIWLLGGLLPFLLLLSGLQRPTEGSEVAIKIDFDFAPGSFDDQYQGCS KQVMEKLTQGDYFTKDIEAQKNYFRMWQKAHLAWLNQGKVL PQNMTTTHAVAILFYTLNSNVHSDFTAM ASVARTPQQYERSFHFHYLTSAILQLLRKDSIMENGTLCYEVHYRTKDVHFNAYTGATIRFGQFLST SLLKEEAQEFNQTLFTIFTCLGAPVQYFSLKKEVLIPPYELFKVINMSYHPRGNWLQLRSTGNLSTYNC QLLKASSKKCIPDPIAIASLSFLT SVIIFSRSRV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_021071
ORF Size:	942 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_021071.2</u>
RefSeq Size:	1386 bp
RefSeq ORF:	945 bp
Locus ID:	420
UniProt ID:	<u>Q93070</u>
Cytogenetics:	12p12.3
MW:	35.9 kDa