

Product datasheet for **RC210366L3V**

SLAMF9 (NM_033438) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLAMF9 (NM_033438) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLAMF9
Synonyms:	CD2F-10; CD2F10; CD84-H1; CD84H1; SF2001
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033438
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210366).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033438.1
RefSeq Size:	1170 bp
RefSeq ORF:	870 bp
Locus ID:	89886
UniProt ID:	Q96A28
Cytogenetics:	1q23.2
Protein Families:	Druggable Genome, Transmembrane
MW:	32.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the signaling lymphocytic activation molecule family. The encoded protein is a cell surface molecule that consists of two extracellular immunoglobulin domains, a transmembrane domain and a short cytoplasmic tail that lacks the signal transduction motifs found in other family members. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Apr 2009]