

Product datasheet for **RC210389L3V**

Dopamine Receptor D1 (DRD1) (NM_000794) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dopamine Receptor D1 (DRD1) (NM_000794) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Dopamine Receptor D1
Synonyms:	DADR; DRD1A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000794
ORF Size:	1338 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210389).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000794.3
RefSeq Size:	3373 bp
RefSeq ORF:	1341 bp
Locus ID:	1812
UniProt ID:	P21728
Cytogenetics:	5q35.2
Domains:	7tm_1
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)

Protein Pathways: Calcium signaling pathway, Gap junction, Neuroactive ligand-receptor interaction

MW: 49.3 kDa

Gene Summary: This gene encodes the D1 subtype of the dopamine receptor. The D1 subtype is the most abundant dopamine receptor in the central nervous system. This G-protein coupled receptor stimulates adenylyl cyclase and activates cyclic AMP-dependent protein kinases. D1 receptors regulate neuronal growth and development, mediate some behavioral responses, and modulate dopamine receptor D2-mediated events. Alternate transcription initiation sites result in two transcript variants of this gene. [provided by RefSeq, Jul 2008]