

## Product datasheet for **RC210395L2V**

### Neurokinin 1 Receptor (TACR1) (NM\_001058) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Neurokinin 1 Receptor (TACR1) (NM_001058) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Neurokinin 1 Receptor
Synonyms:	NK1R; NKIR; SPR; TAC1R
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001058
ORF Size:	1221 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210395).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001058.2</a>
RefSeq Size:	4795 bp
RefSeq ORF:	1224 bp
Locus ID:	6869
UniProt ID:	<a href="#">P25103</a>
Cytogenetics:	2p12
Domains:	7tm_1
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)

**Protein Pathways:** Calcium signaling pathway, Neuroactive ligand-receptor interaction

**MW:** 46.3 kDa

**Gene Summary:** This gene belongs to a gene family of tachykinin receptors. These tachykinin receptors are characterized by interactions with G proteins and contain seven hydrophobic transmembrane regions. This gene encodes the receptor for the tachykinin substance P, also referred to as neurokinin 1. The encoded protein is also involved in the mediation of phosphatidylinositol metabolism of substance P. [provided by RefSeq, Sep 2008]