

## Product datasheet for **RC210396L3V**

### HOXB2 (NM\_002145) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HOXB2 (NM_002145) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HOXB2
Synonyms:	Hox-2.8; HOX2; HOX2H; K8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002145
ORF Size:	1068 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210396).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002145.3</a> , <a href="#">NP_002136.1</a>
RefSeq Size:	1614 bp
RefSeq ORF:	1071 bp
Locus ID:	3212
UniProt ID:	<a href="#">P14652</a>
Cytogenetics:	17q21.32
MW:	37.9 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the Antp homeobox family and encodes a nuclear protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development. Increased expression of this gene is associated with pancreatic cancer. [provided by RefSeq, Jul 2008]