

Product datasheet for **RC210404L3V**

Proteinase Activated Receptor 4 (F2RL3) (NM_003950) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Proteinase Activated Receptor 4 (F2RL3) (NM_003950) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Proteinase Activated Receptor 4
Synonyms:	PAR4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003950
ORF Size:	1155 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210404).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003950.2
RefSeq Size:	2731 bp
RefSeq ORF:	1158 bp
Locus ID:	9002
UniProt ID:	Q96RI0
Cytogenetics:	19p13.11
Domains:	7tm_1



[View online »](#)

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Neuroactive ligand-receptor interaction

MW: 41.13 kDa

Gene Summary: This gene encodes a member of the protease-activated receptor subfamily, part of the G-protein coupled receptor 1 family of proteins. The encoded receptor is proteolytically processed to reveal an extracellular N-terminal tethered ligand that binds to and activates the receptor. This receptor plays a role in blood coagulation, inflammation and response to pain. Hypomethylation at this gene may be associated with lung cancer in human patients. [provided by RefSeq, Sep 2016]