

## Product datasheet for **RC210417L3V**

### TRAF3 (NM\_145725) Human Tagged ORF Clone Lentiviral Particle

#### Product data:


Product Type:	Lentiviral Particles
Product Name:	TRAF3 (NM_145725) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRAF3
Synonyms:	CAP-1; CAP1; CD40bp; CRAF1; IIAE5; LAP1; RNF118
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145725
ORF Size:	1704 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210417).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145725.1</a>
RefSeq Size:	7793 bp
RefSeq ORF:	1707 bp
Locus ID:	7187
UniProt ID:	<a href="#">Q13114</a>
Cytogenetics:	14q32.32
Protein Families:	Druggable Genome



[View online »](#)

<b>Protein Pathways:</b>	Pathways in cancer, RIG-I-like receptor signaling pathway, Small cell lung cancer, Toll-like receptor signaling pathway
<b>MW:</b>	64.5 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the TNF receptor associated factor (TRAF) protein family. TRAF proteins associate with, and mediate the signal transduction from, members of the TNF receptor (TNFR) superfamily. This protein participates in the signal transduction of CD40, a TNFR family member important for the activation of the immune response. This protein is found to be a critical component of the lymphotoxin-beta receptor (LTbetaR) signaling complex, which induces NF-kappaB activation and cell death initiated by LTbeta ligation. Epstein-Barr virus encoded latent infection membrane protein-1 (LMP1) can interact with this and several other members of the TRAF family, which may be essential for the oncogenic effects of LMP1. The protein also plays a role in the regulation of antiviral response. Mutations in this are associated with Encephalopathy, acute, infection-induced, herpes-specific 5. [provided by RefSeq, Jul 2020]

### Product images:



[RC210417L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC210417L3V particle to overexpress human TRAF3-Myc-DDK fusion protein.