

Product datasheet for **RC210446L3V**

KCNK10 (NM_138318) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KCNK10 (NM_138318) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KCNK10
Synonyms:	K2p10.1; PPP1R97; TREK-2; TREK2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138318
ORF Size:	1629 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210446).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138318.1
RefSeq Size:	7198 bp
RefSeq ORF:	1632 bp
Locus ID:	54207
UniProt ID:	P57789
Cytogenetics:	14q31.3
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
MW:	60.1 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the family of potassium channel proteins containing two pore-forming P domains. This channel is an open rectifier which primarily passes outward current under physiological K⁺ concentrations, and is stimulated strongly by arachidonic acid and to a lesser degree by membrane stretching, intracellular acidification, and general anaesthetics. Several alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Sep 2008]