

Product datasheet for **RC210467L4V**

LBH (NM_030915) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LBH (NM_030915) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LBH
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030915
ORF Size:	315 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210467).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030915.1
RefSeq Size:	2956 bp
RefSeq ORF:	318 bp
Locus ID:	81606
UniProt ID:	Q53QV2
Cytogenetics:	2p23.1
MW:	12.2 kDa
Gene Summary:	Transcriptional activator which may act in mitogen-activated protein kinase signaling pathway.[UniProtKB/Swiss-Prot Function]



[View online »](#)