

Product datasheet for **RC210468L4V**

ANKRD2 (NM_020349) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ANKRD2 (NM_020349) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ANKRD2
Synonyms:	ARPP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020349
ORF Size:	1080 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210468).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020349.2
RefSeq Size:	1520 bp
RefSeq ORF:	1083 bp
Locus ID:	26287
UniProt ID:	Q9GZV1
Cytogenetics:	10q24.2
MW:	39.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that belongs to the muscle ankyrin repeat protein (MARP) family. A similar gene in rodents is a component of a muscle stress response pathway and plays a role in the stretch-response associated with slow muscle function. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2014]