

## Product datasheet for **RC210470L1V**

### **PVRIG (NM\_024070) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**


Product Type:	Lentiviral Particles
Product Name:	PVRIG (NM_024070) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PVRIG
Synonyms:	C7orf15; CD112R
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_024070
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210470).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024070.3</a> , <a href="#">NP_076975.2</a>
RefSeq Size:	1587 bp
RefSeq ORF:	981 bp



[View online »](#)

Locus ID:	79037
UniProt ID:	<a href="#">Q6DKI7</a>
Cytogenetics:	7q22.1
Protein Families:	Transmembrane
MW:	34.3 kDa
Gene Summary:	Cell surface receptor for NECTIN2. May act as a coinhibitory receptor that suppresses T-cell receptor-mediated signals. Following interaction with NECTIN2, inhibits T-cell proliferation. Competes with CD226 for NECTIN2-binding.[UniProtKB/Swiss-Prot Function]

### Product images:



[RC210470L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC210470L1V particle to overexpress human PVRIG-Myc-DDK fusion protein.