

Product datasheet for RC210470L4

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

PVRIG (NM_024070) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: PVRIG (NM_024070) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: PVRIG

Synonyms: C7orf15; CD112R

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC210470).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_024070

ORF Size: 978 bp



OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at customport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

Cytogenetics:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM 024070.3, NP 076975.2

7q22.1

 RefSeq Size:
 1587 bp

 RefSeq ORF:
 981 bp

 Locus ID:
 79037

 UniProt ID:
 Q6DKI7

Protein Families: Transmembrane

MW: 34.3 kDa


Gene Summary: Cell surface receptor for NECTIN2. May act as a coinhibitory receptor that suppresses T-cell

receptor-mediated signals. Following interaction with NECTIN2, inhibits T-cell proliferation.


Competes with CD226 for NECTIN2-binding.[UniProtKB/Swiss-Prot Function]



Product images:



Circular map for RC210470L4



Double digestion of RC210470L4 using Sgfl and Mlul $\,$