

Product datasheet for **RC210472L4V**

CFHL1 (CFHR1) (NM_002113) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CFHL1 (CFHR1) (NM_002113) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CFHL1
Synonyms:	CFHL; CFHL1; CFHL1P; CFHR1P; FHR-1; FHR1; H36; H36-1; H36-2; HFL1; HFL2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002113
ORF Size:	990 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210472).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002113.2
RefSeq Size:	1312 bp
RefSeq ORF:	993 bp
Locus ID:	3078
UniProt ID:	Q03591
Cytogenetics:	1q31.3
Domains:	CCP
Protein Families:	Secreted Protein



[View online »](#)

MW: 37.7 kDa

Gene Summary: This gene encodes a secreted protein belonging to the complement factor H protein family. It binds to *Pseudomonas aeruginosa* elongation factor Tuf together with plasminogen, which is proteolytically activated. It is proposed that Tuf acts as a virulence factor by acquiring host proteins to the pathogen surface, controlling complement, and facilitating tissue invasion. Mutations in this gene are associated with an increased risk of atypical hemolytic-uremic syndrome. [provided by RefSeq, Oct 2009]