

Product datasheet for **RC210475**

SELS (SELENOS) (NM_203472) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SELS (SELENOS) (NM_203472) Human Tagged ORF Clone
Symbol: SELS
Synonyms: AD-015; ADO15; SBB18; SELS; SEPS1; VIMP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC210475 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAACGCCAAGAGGAGTCTCTGTCCGCGCGCCGGCCCTGGAGACCGAGGGGCTGCGCTTCTGCACA
 CCACGGTGGGCTCCCTGCTGGCCACCTATGGCTGGTACATCGTCTTCAGCTGCATCCTTCTCTACGTGGT
 CTTTCAGAAGCTTCCGCCCGGCTAAGAGCCTTGAGGCAGAGGCAGCTGGACCGAGCTGCGGCTGCTGTG
 GAACCTGATGTTGTTGTTAAACGACAAGAAGCTTTAGCAGCTGCTCGACTGAAAAATGCAAGAAGAATAA
 ATGCGCAAGTTGAAAAGCATAAGGAAAACTGAAACAATTGAAGAAGAAAAAGGAGACAGAGAATTGA
 AATGTGGGACAGCATGCAAGAAGGAAAAAGTTACAAAGGAAATGCAAAGAAGCCCGAGGAGGAGCAGT
 CCTGGGCTTCCACTTCTGTCTGCTGAAACGGAATCGGACAGAAAGCCTTTGCGGGGAGGAGTTATA
 ACCCGTTGTCTGGTGAAGGAGGCGGAGCTTGTCTGGAGACCTGGACGAGAGGCCCGTCATCTGGCGG
 ATGAGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210475 protein sequence
 Red=Cloning site Green=Tags(s)

MERQEESLSARPALETEGLRFLHTTVGSLLATYGWYIVFSCILLYVVFQKLSARLRALRQRQLDRAAAAV
 EPDVVVKRQEALAAARLKMQEELNAQVEKHKEKQLEEEKRRQK IEMWDSMQEGKSYKGNAKKPQEEDS
 PGPSTSSVLKRKSDRKLPRGGYNPLSGEGGGACSWRPGRRGPSSGG*G

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6547_a01.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_203472

ORF Size: 567 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone encodes a selenoprotein containing the rare amino acid selenocysteine (Sec). Sec is encoded by UGA codon, which normally signals translational termination. Expression of this clone is not guaranteed due to the nature of selenoproteins.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_203472.1](#), [NP_982298.1](#)

RefSeq Size: 1417 bp

RefSeq ORF: 564 bp

Locus ID: 55829


UniProt ID: [Q9BQE4](#)

Cytogenetics: 15q26.3

Protein Families: Druggable Genome

Gene Summary: This gene encodes a transmembrane protein that is localized in the endoplasmic reticulum (ER). It is involved in the degradation process of misfolded proteins in the ER, and may also have a role in inflammation control. This protein is a selenoprotein, containing the rare amino acid selenocysteine (Sec). Sec is encoded by the UGA codon, which normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, designated the Sec insertion sequence (SECIS) element, that is necessary for the recognition of UGA as a Sec codon, rather than as a stop signal. Two additional phylogenetically conserved stem-loop structures (Stem-loop 1 and Stem-loop 2) in the 3' UTR of this mRNA have been shown to function as modulators of Sec insertion. An alternatively spliced transcript variant, lacking the SECIS element and encoding a non-Sec containing shorter isoform, has been described for this gene (PMID:23614019). [provided by RefSeq, Jul 2017]

Product images:



Circular map for RC210475