

## Product datasheet for **RC210482L4V**

### NDUFA4 (NM\_002489) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NDUFA4 (NM_002489) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NDUFA4
Synonyms:	CI-9k; CI-MLRQ; COXFA4; MC4DN21; MISTR1; MLRQ; MRCAF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002489
ORF Size:	243 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210482).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002489.2</a>
RefSeq Size:	2058 bp
RefSeq ORF:	246 bp
Locus ID:	4697
UniProt ID:	<a href="#">O00483</a>
Cytogenetics:	7p21.3
Protein Families:	Transmembrane



[View online »](#)

<b>Protein Pathways:</b>	Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
<b>MW:</b>	9.4 kDa
<b>Gene Summary:</b>	The protein encoded by this gene belongs to the complex I 9kDa subunit family. Mammalian complex I of mitochondrial respiratory chain is composed of 45 different subunits. This protein has NADH dehydrogenase activity and oxidoreductase activity. It transfers electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone. [provided by RefSeq, Jul 2008]