

Product datasheet for **RC210489A1V**

Human Bcl 7A (BCL7A) (NM_001024808) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Bcl 7A (BCL7A) (NM_001024808) AAV Particle
Tag: Myc-DDK
Symbol: Bcl 7A
Synonyms: BCL7
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210489 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCGGGCAGGTCGGTTCGAGCCGAGACGAGGAGCCGGGCCAAAGATGATATCAAGAGGGTCATGGCGG
 CGATCGAGAAAGTGCCTAAATGGGAGAAGAAATGGGTGACCGTTGGTGACACATCCCTACGAATCTACAA
 ATGGGTCCCTGTGACGGAGCCCAAGGTTGATGACAAAAACAAGAATAAGAAAAAGGCAAGGACGAGAAG
 TGTGGCTCAGAGGTGACCACTCCGGAGAACAGTTCTCCCCAGGGATGATGGACATGCATGACGATAACA
 GCAACCCAGAGCTCCATCGCAGATGCCTCCCCATCAACAGGAGAACAGCAGCAACTCCAGCCCCGCTCC
 AGAGCCCAACTCGGCTGTGCCAGCGACGGCACCGAGGCCAAGGTGGATGAGGCCAGGCTGATGGGAAG
 GAGCACCCAGGAGCTGAAGATGCTTCTGATGAGCAGAATTCACAGTCCTCGATGGAACATTTCATGAACA
 GCTCAGAGAAAGTAGATCGGCAGCCGTCTGGAGACTCGGGTCTGGCCGAGAGACGTCTGCAATCTCTCA
 GGATTTGGAAGGAGTGCCACCCTCTAAAAAGATGAACTGGAGGCCTCTCAACAAAACCCGAAGAGATG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210489 protein sequence
 Red=Cloning site Green=Tags(s)

MSGRSVRAETRSRAKDDIKRVMAAIEKVRKWEKKWVTVGDTSLRIYKWPVTEPKVDDKNKNKKKGKDEK
 CGSEVTPPENSSSPGMMDMHDDNSNQSSIADASPIKQENSNSSPAPEPNSAVPSDGTEAKVDEAQADGK
 EHPGAEASDEQNSQSSMEHSMNSSEKVDQRQPSGDSGLAAETSAISQDLEGVPPSKMKLEASQQNSEEM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001024808
ORF Size:	630 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001024808.1</u>
RefSeq Size:	3651 bp
RefSeq ORF:	633 bp
Locus ID:	605
UniProt ID:	<u>Q4VC05</u>
Cytogenetics:	12q24.31
MW:	22.8 kDa