

Product datasheet for **RC210495L3V**

SLC2A11 (NM_030807) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC2A11 (NM_030807) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC2A11
Synonyms:	GLUT10; GLUT11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030807
ORF Size:	1509 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210495).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030807.2
RefSeq Size:	3133 bp
RefSeq ORF:	1512 bp
Locus ID:	66035
UniProt ID:	Q9BYW1
Cytogenetics:	22q11.23
Domains:	sugar_tr
Protein Families:	Transmembrane



[View online »](#)

MW: 54.4 kDa

Gene Summary: This gene belongs to a family of proteins that mediate the transport of sugars across the cell membrane. The encoded protein transports glucose and fructose. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]