

Product datasheet for **RC210529L3V**

CTBP1 (NM_001328) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CTBP1 (NM_001328) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CTBP1
Synonyms:	BARS; HADDTS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001328
ORF Size:	1320 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210529).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001328.2
RefSeq Size:	2288 bp
RefSeq ORF:	1323 bp
Locus ID:	1487
UniProt ID:	Q13363
Cytogenetics:	4p16.3
Domains:	2-Hacid_DH, 2-Hacid_DH_C



[View online »](#)

Protein Pathways:	Chronic myeloid leukemia, Notch signaling pathway, Pathways in cancer, Wnt signaling pathway
MW:	47.5 kDa
Gene Summary:	This gene encodes a protein that binds to the C-terminus of adenovirus E1A proteins. This phosphoprotein is a transcriptional repressor and may play a role during cellular proliferation. This protein and the product of a second closely related gene, CTBP2, can dimerize. Both proteins can also interact with a polycomb group protein complex which participates in regulation of gene expression during development. Alternative splicing of transcripts from this gene results in multiple transcript variants. [provided by RefSeq, Jul 2008]