

Product datasheet for **RC210537L3V**

FAM40A (STRIP1) (NM_033088) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM40A (STRIP1) (NM_033088) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM40A
Synonyms:	FAM40A; FAR11A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033088
ORF Size:	2511 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210537).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033088.2
RefSeq Size:	3291 bp
RefSeq ORF:	2514 bp
Locus ID:	85369
UniProt ID:	Q5VSL9
Cytogenetics:	1p13.3
MW:	95.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the striatin-interacting phosphatase and kinase complex, which is involved in localization of the Golgi body. The encoded protein participates in cytoskeletal organization. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2012]