

## Product datasheet for **RC210545L3V**

### **ANO6 (NM\_001025356) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	ANO6 (NM_001025356) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ANO6
Synonyms:	BDPLT7; SCTS; TMEM16F
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025356
ORF Size:	2730 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210545).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001025356.2</a>
RefSeq Size:	6133 bp
RefSeq ORF:	2733 bp
Locus ID:	196527
UniProt ID:	<a href="#">Q4KMQ2</a>
Cytogenetics:	12q12
Protein Families:	Transmembrane
MW:	106.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a multi-pass transmembrane protein that belongs to the anoctamin family. This protein is an essential component for the calcium-dependent exposure of phosphatidylserine on the cell surface. The scrambling of phospholipid occurs in various biological systems, such as when blood platelets are activated, they expose phosphatidylserine to trigger the clotting system. Mutations in this gene are associated with Scott syndrome. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Mar 2011]