

## Product datasheet for **RC210565L3V**

### **KLF10 (NM\_005655) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KLF10 (NM_005655) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLF10
Synonyms:	EGR-alpha; EGRA; TIEG; TIEG1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005655
ORF Size:	1440 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210565).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005655.1</a>
RefSeq Size:	3115 bp
RefSeq ORF:	1443 bp
Locus ID:	7071
UniProt ID:	<a href="#">Q13118</a>
Cytogenetics:	8q22.3
Domains:	zf-C2H2
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 52.6 kDa

**Gene Summary:** This gene encodes a member of a family of proteins that feature C2H2-type zinc finger domains. The encoded protein is a transcriptional repressor that acts as an effector of transforming growth factor beta signaling. Activity of this protein may inhibit the growth of cancers, particularly pancreatic cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2013]