

## Product datasheet for **RC210568L4V**

### DMT1 (SLC11A2) (NM\_000617) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DMT1 (SLC11A2) (NM_000617) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DMT1
Synonyms:	AHMIO1; DCT1; DMT1; NRAMP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000617
ORF Size:	1683 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210568).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000617.1</a>
RefSeq Size:	4164 bp
RefSeq ORF:	1686 bp
Locus ID:	4891
UniProt ID:	<a href="#">P49281</a>
Cytogenetics:	12q13.12
Domains:	Nramp
Protein Families:	Transmembrane



[View online »](#)

**Protein Pathways:** Lysosome

**MW:** 61.5 kDa

**Gene Summary:** This gene encodes a member of the solute carrier family 11 protein family. The product of this gene transports divalent metals and is involved in iron absorption. Mutations in this gene are associated with hypochromic microcytic anemia with iron overload. A related solute carrier family 11 protein gene is located on chromosome 2. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Apr 2010]