

Product datasheet for **RC210587A1V**

Human PLSCR3 (NM_020360) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PLSCR3 (NM_020360) AAV Particle
Tag: Myc-DDK
Symbol: PLSCR3
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210587 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGGCTACTTGCCCCCAAAGGCTACGCCCTTCGCCCCACCTCCCTACCCTGTCACCCCTGGGT
ACCCGGAGCCGGCGCTACATCCTGGGCCCGGGCAGGCGCCAGTGCCCGCCAGGTACCTGCCCCAGCTCC
CGGCTTCGCCCTCTCCCTCGCCTGGCCCCGTGGCCTGGGGTCTGCTGCCCTTCTTGCCACTGCCA
GGGTGCCTTCTGGCCTCGAATTCCTGGTGCAGATTGATCAGATTTTGATTACCAGAAGGCTGAGCGAG
TGAAACGTTCTAGGCTGGGAGACCTGTAATCGGTATGAACTGCGCTCTGGGGCCGGGCAGCCCTGGG
TCAGGCGGCCGAGGAGAGCAACTGCTGCGCCGCTGTGCTGTGGCGCCCGCCGCGCTGCGTGTCCGC
CTGGCCGACCCCGGGGACCGTGAGGTGCTGCGTTTGCTCCGCCGCTGCACTGTGGCTGCAGCTGCTGCC
CCTGTGGCCTCCAGGAGATGGAAGTACAGGCTCCACCAGGCACCACTTGCCACGTGCTACAGACCTG
GCATCCCTTCTCCCAAGTTCTCCATCCAGGATGCCGATCGCCAGACAGTCTTGCGAGTGGTGGGGCC
TGCTGGACCTGTGGCTGTGGCACAGACCAACTTTGAGGTGAAGACTCGGGATGAATCCCGCAGTGTGG
GCCGATCAGCAAGCAGTGGGGGGCCCTGGTCCGAGAAGCCCTCACAGATGCAGATGACTTTGGCCTACA
GTTCCCGCTGGACCTGGATGTGAGGGTGAAGGCTGTGCTGCTGGGAGCCACATTCCTCATTGACTACATG
TTCTTTGAGAAGCGAGGAGGCGCTGGGCCCTCTGCCGTCACCAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210587 protein sequence Red=Cloning site Green=Tags(s) MAGYLPPKGYAPSPPPYPVTPGYPEPALHPGPGQAPVPAQVPAPAPGFALFPSPGPVALGSAAPFLPLP GVP SGLEFLVQIDQILIHQKAERVETFLGWETCNRYELRSGAGQPLGQAAEESNCCARLCCGARRPLRVR LADPGDREVLRLRLHCGCSCCPCGLQEMEVQAPP GTTIGHVLQTHWPF LPKFSIQDADRQTVLRVVG P CWT CGCGTDTNFEVKTRDESRSVGRISKQWGGLVREALTDADDFGLQFPLDL D VRVKAVLLGATFLIDYM FFEKRGGAGPSAVTS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_020360
ORF Size:	885 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_020360.2
RefSeq Size:	2089 bp
RefSeq ORF:	888 bp
Locus ID:	57048
UniProt ID:	Q9NRY6
Cytogenetics:	17p13.1
MW:	31.6 kDa