

Product datasheet for **RC210595L3V**

STEAP3 (NM_001008410) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	STEAP3 (NM_001008410) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STEAP3
Synonyms:	AHMIO2; dudlin-2; dudulin-2; pHyde; STMP3; TSAP6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001008410
ORF Size:	1464 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210595).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008410.1 , NP_001008410.1
RefSeq Size:	3870 bp
RefSeq ORF:	1467 bp
Locus ID:	55240
UniProt ID:	Q658P3
Cytogenetics:	2q14.2
Protein Families:	Transmembrane
Protein Pathways:	p53 signaling pathway



[View online »](#)

MW: 54.6 kDa

Gene Summary: This gene encodes a multipass membrane protein that functions as an iron transporter. The encoded protein can reduce both iron (Fe³⁺) and copper (Cu²⁺) cations. This protein may mediate downstream responses to p53, including promoting apoptosis. Deficiency in this gene can cause anemia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2015]