

## Product datasheet for RC210606L3

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## alpha Synuclein (SNCA) (NM\_000345) Human Tagged Lenti ORF Clone

### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** alpha Synuclein (SNCA) (NM\_000345) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: alpha Synuclein

Synonyms: NACP; PARK1; PARK4; PD1

Mammalian Cell

Selection:

Puromycin

**Vector:** pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RC210606).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_000345

ORF Size: 420 bp





#### alpha Synuclein (SNCA) (NM\_000345) Human Tagged Lenti ORF Clone - RC210606L3

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 000345.2</u>

RefSeq Size: 1543 bp RefSeq ORF: 423 bp

Locus ID: 6622

UniProt ID: P37840

Cytogenetics: 4q22.1

**Domains:** Synuclein

**Protein Families:** Druggable Genome

**Protein Pathways:** Alzheimer's disease, Parkinson's disease

MW: 14.3 kDa

**Gene Summary:** Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-

synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively. SNCA may serve to integrate presynaptic signaling and


membrane trafficking. Defects in SNCA have been implicated in the pathogenesis of

Parkinson disease. SNCA peptides are a major component of amyloid plaques in the brains of patients with Alzheimer's disease. Alternatively spliced transcripts encoding different isoforms


have been identified for this gene. [provided by RefSeq, Feb 2016]



# **Product images:**



Circular map for RC210606L3



Double digestion of RC210606L3 using Sgfl and Mlul