

Product datasheet for **RC210611A1V**

Human SNX11 (NM_013323) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNX11 (NM_013323) AAV Particle
Tag: Myc-DDK
Symbol: SNX11
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210611 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGCTTTTGGTGTAGGATGTCGGAGAACCAAGAACAGGAGGAGGTGATTACAGTGCCTGTTTCAGGACC
CCCGAGTGCAGAATGAGGGCTCCTGGAACCTTATGTGGATTATAAGATATTCCTCCATACCAACAGCAA
AGCCTTTACTGCCAAGACTTCCTGTGTGCGGCGCCGCTACCGTGAGTTCGTGTGGCTGAGAAAGCAGCTA
CAGAGAAATGCTGGTTTGGTGCCTGTTCTGAACCTCCTGGGAAGTCAACCTTCTCGGCACCTCAGATG
AGTTCATTGAGAAGCGACGACAAGGTCTGCAGCACTTCCTTGAAAAGGTCTGCAGAGTGTGGTTCTCCT
GTCAGACAGCCAGTTGCACCTATTCCTGCAAAGCCAGCTCTCGGTGCCTGAGATAGAAGCCTGTGTCCAG
GGCCGAAGTACCATGACTGTGTCTGATGCCATTCTTCGATATGCTATGTCAAACCTGTGGCTGGGCCCAGG
AAGAGAGGCAGAGCTCTTCTCACCTGGCTAAAGGAGACCAGCCTAAGAGTTGCTGCTTTCTTCCAAGATC
GGGTAGGAGGAGCTCTCCCTCACCGCCTCCCAGTGAAGAAAAGGACCATTTAGAAGTGTGGGCTCCAGTT
GTTGACTCTGAGGTTCTTCTTGAAAAGTCCCACTCTCCACCCCTCTCCTCACCATTATGCTGTGATT
TTGGAAGACCCAAAGAGGGAACCTCCACTCTTCAGTCTGTGAGGAGGGCTGTGGGAGGAGATCATGCTGT
GCCTTTGGACCCTGGTCAGTTAGAAACAGTTTTGGAAAAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210611 protein sequence Red=Cloning site Green=Tags(s) MGFWCRMSENQEEVITVRVQDPRVQNEGSWNSYVDYKIFLHTNSKAFTAKTSCVRRRYREFVWLRKQL QRNAGLVPVPELPGKSTFFGTSDEFIEKRRQGLQHFLEKVLQSVVLLSDSQLHLFLQSQLSVPEIEACVQ GRSTMTVSDAILRYAMSNCGWAQEERQSSSHLAKGDQPKSCCFLPRSGRRSSPSPPPSEEKDHLEWAPV VDSEVPSLESPTLPPLSSPLCCDFGRPKEGTSTLQSVRRAVGGDHAVPLDPGQLETVLEK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_013323
ORF Size:	810 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_013323.2
RefSeq Size:	2315 bp
RefSeq ORF:	813 bp
Locus ID:	29916
UniProt ID:	Q9Y5W9
Cytogenetics:	17q21.32
MW:	30.4 kDa