

## Product datasheet for **RC210621L1V**

### TEAD3 (NM\_003214) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TEAD3 (NM_003214) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TEAD3
Synonyms:	DTEF-1; ETEFR-1; TEAD-3; TEAD5; TEF-5; TEF5
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003214
ORF Size:	1305 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210621).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003214.2</a>
RefSeq Size:	3023 bp
RefSeq ORF:	1308 bp
Locus ID:	7005
UniProt ID:	<a href="#">Q99594</a>
Cytogenetics:	6p21.31
Domains:	TEA
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 48.7 kDa

**Gene Summary:** This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and is involved in the transactivation of the chorionic somatomammotropin-B gene enhancer. Translation of this protein is initiated at a non-AUG (AUA) start codon. [provided by RefSeq, Jul 2008]