

Product datasheet for **RC210664A1V**

Human Serum Amyloid A (SAA1) (NM_000331) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Serum Amyloid A (SAA1) (NM_000331) AAV Particle
Tag: Myc-DDK
Symbol: Serum Amyloid A
Synonyms: PIG4; SAA; SAA2; TP53I4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210664 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAGCTTCTCACGGGCTGGTTTTCTGCTCCTTGGTCTGGGTGTCAGCAGCCGAAGCTTCTTTTCGT
 TCCTTGGCGAGGCTTTTGATGGGGCTCGGGACATGTGGAGAGCCTACTCTGACATGAGAGAAGCCAATTA
 CATCGGCTCAGACAAATACTTCCATGCTCGGGGAACATGATGCTGCCAAAAGGGGACCTGGGGGTGTC
 TGGGCTCAGAAAGCGATCAGCGATGCCAGAGAGAATATCCAGAGATTCTTTGGCCATGGTGCGGAGGACT
 CGCTGGCTGATCAGGCTGCCAATGAATGGGGCAGGAGTGGCAAAGACCCCAATCACTTCCGACCTGCTGG
 CCTGCCTGAGAAATAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210664 protein sequence
 Red=Cloning site Green=Tags(s)

MKLLTGLVFCSLVLGVSSRSFFSFLGEAFDGDARDMWRAYSDMREANYIGSDKYFHARGNYDAAKRGPGGV
 WAAEAISDARENIQRFFGHGAEDSLADQAANEWGRSGKDPNHFRPAGLPEKY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_000331



[View online »](#)

ORF Size:	366 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_000331.2</u>
RefSeq Size:	678 bp
RefSeq ORF:	369 bp
Locus ID:	6288
UniProt ID:	<u>P02735</u>
Cytogenetics:	11p15.1
MW:	13.5 kDa