

## Product datasheet for **RC210670**

### **TSNAXIP1 (NM\_018430) Human Tagged ORF Clone**

#### **Product data:**

|                           |                                             |
|---------------------------|---------------------------------------------|
| Product Type:             | Expression Plasmids                         |
| Product Name:             | TSNAXIP1 (NM_018430) Human Tagged ORF Clone |
| Tag:                      | Myc-DDK                                     |
| Symbol:                   | TSNAXIP1                                    |
| Synonyms:                 | TXI1                                        |
| Mammalian Cell Selection: | Neomycin                                    |
| Vector:                   | pCMV6-Entry (PS100001)                      |
| E. coli Selection:        | Kanamycin (25 ug/mL)                        |



[View online »](#)

**ORF Nucleotide  
Sequence:**

&gt;RC210670 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGC**C

ATGGGTGGGCACCTGTCCCCATGGCCACATACACCAGTGGCCAGACCATTTTGCAAATCGAAAACCTT  
 GTTCAGATGACTACCGGAAGCGAGTAGGGAGCTGCCAGCAGACCCCTTTCGCACTGCCAAGCCCCAGTA  
 CTTGGAGGAAGTGGAAAACTACCTACGCAAGGAGCTCCTCCTGCTGGACCTGGGCACAGATTCCACCCAG  
 GAACTAAGGCTGCAGCCTTACAGAGAGATCTTTGAGTTCTTCATAGAGGACTTCAAAACGTACAAGCCAT  
 TACTATCTCCATCAAGAATGCGTATGAGGGGATGCTGGCCACCAAAGGAGAAGATTCGGGCTCTGGA  
 GCCCTGAAGGCCAAGCTTGTCACTGTGAATGAGGACTGCAATGAGAGGATCCTGGCCATGAGAGCTGAG  
 GAGAAATATGAAATCTCCTGCTCAAGAAAGAGAAGATGAATTTGCTAAACTCATCGACAAAAAGAATG  
 AGGAGAAGATTTTCATTGCAGAGCGAGGTGACCAAACTGAGGAAGAAGCTTGGCTGAGGAGTACCTGCACTA  
 CCTCAGTGAGCGAGATGCCTGTAAAGATCCTCATCGCAGACCTGAATGAGCTGCGGTACCAGCGGGAGGAC  
 ATGTCATTAGCCAGTCGCCAGGCATCTGGGGGAGGACCCTGTGAAGTTAACCTGGCTCTTAAGATGA  
 CCCGGCAAGACCTGACCCGCACGCAGATGGAACCTCAACAACATGAAGGCCAACTTTGGAGATGTGGTCCC  
 CAGGAGGGACTTTGAAATGCAGGAGAAGACCAACAAGGATCTTCAGGAGCAGCTGGACACCTGAGAGCC  
 AGCTACGAGGAGGTTTCGAAGGAGCATGAGATCCTCATGCAGCTGCACATGAGCAGCTGAAGGAACGGG  
 ACCAATTTCTTCTGAGCTGCAGGAGATCCAGCGCACTTCCACGCCGCGGCCTGACTGGACCAAGTGCAA  
 AGATGTGGTGGCTGGGGGCCAGAGCGCTGGCAGATGCTGGCTGAGGGCAAGAACAGCGACCAGCTGGTG  
 GACGTGCTCCTGGAAGAGATTGGTTCGGGGCTGCTGCGGGAGAAAGACTTCTTCCCTGGTCTGGGCTATG  
 GGGAAGCCATCCCTGCTTTTCTCGGTTTGATGGCCTCGTGGAGAACAGAAGCCAAGCAAGAAGGACGT  
 GGTCAACCTCCTCAAGGATGCCTGGAAGGAACGTCTTGCTGAGGAGCAGAAAGAGACGTTCCAGATTTT  
 TTCTTCAATTTCTGAGCATCGCTTTGGGCCAGTGATGCCATGGCCTGGGCTTATACTATTTTGA  
 ATATCAAGATCTTCACTCCAACGAGGTTATGAGTCAGTTCTATGCAGTCTTGATGGGAAAGCGGAGTGA  
 GAATGTGTATGTCACCCAGAAGGAGACAGTAGCCAGCTGCTGAAGGAGATGACAAATGCTGACAGTCAG  
 AACGAGGGGCTACTAACCATGGAGCAGTTCAACACTGTCCTCAAGAGTACCTTCCCTCTCAAGACAGAAG  
 AGCAAATCCAGGAGCTGATGGAGGCAGGGGCTGGCATCCAGCAGCAGCAATGCAGACTTGCTCAACTA  
 CCGCTCACTGTTTATGGAGGATGAGGAGGGCCAGAGTGAGCCCTTGTGCAAAAACCTCTGGGAACAATAC  
 ATGGATGAGAAGGACGAGTACTTACAGCAGCTAAAGCAGGAGCTTGGCATAGAACTCCATGAGGAAGTGA  
 CTCTGCCCAAGCTGCGAGGGGGCTGATGACCATCGACCCAGCCTGGACAAGCAGACAGTGAACACCTA  
 CATGAGCCAGGCTTCCAGCTCCCTGAGTCGGAATGCCAGAGGAGGGTACGAGAAGGAAGAAGCCGTG  
 GTGGAAATCCTCCAGACTGCCCTGGAGCGGCTTCAGGTGATTGACATCAGGCGTGTGGGACCTCGAGAGC  
 CAGAGCCTGCAAGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC210670 protein sequence  
 Red=Cloning site Green=Tags(s)


MGGHLSWPPTYTSGQTILQNRKPCSDDYRKRVGSCQHPFRTAKPQYLEELENYLRKELLLDLGTDSTQ  
 ELRLQPYREIFEFFIEDFKTYKPLLSSIKNAYEGMLAHQREKIRALEPLKAKLVTVNEDCNERILAMRAE  
 EKYEISLLKKEKMNLKLDKKNEEKISLQSEVTKLRKNLAEEYLHYLSERDACKILIALDLNELRYQRED  
 MSLAQSPGIWGEDPVKLTALKMTRQDLTRTQMELNNMKANFGDVVPRDFEMQEKTNKDLQEQLDTLRA  
 SYEEVRKEHEILMQLHMSTLKERDQFFSELQEIQRTSTPRPDWTKCKDVVAGGPERWQMLAEGKNSDQLV  
 DVLLLEEIGSGLLREKDFPGLGYGEAIPAFLRFDGLVENKKPSKKDVVNLKDAWKERLAEEQKETFPDF  
 FNFLEHRFGPSDAMAWAYTIFENIKIFHSNEVMSQFYAVLMGKRSENVYVTQKETVAQLLKEMTNADSQ  
 NEGLLTMEQFNTVLKSTFPLKTEEQIQELMEAGGWHPPSSNADLLNYSRLFMEDEEGQSEPFVQKLWEQY  
 MDEKDEYLQQLKQELGIELHEEVTLPKLRGGLMTIDPSLDKQTVNTYMSQAFQLPESEMPEEGDEKEEAV  
 VEILQTALERLQVIDIRRVGPPEPEPAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6585\\_e10.zip](https://cdn.origene.com/chromatograms/mk6585_e10.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_018430

**ORF Size:** 1974 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_018430.4](#)

**RefSeq Size:** 2692 bp

**RefSeq ORF:** 1977 bp

**Locus ID:** 55815


**UniProt ID:** [Q2TAA8](#)

**Cytogenetics:** 16q22.1


**MW:** 76.8 kDa

**Gene Summary:** Possible role in spermatogenesis.[UniProtKB/Swiss-Prot Function]


## Product images:



Circular map for RC210670



Western blot validation of overexpression lysate (Cat# [LY413059]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC210670 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified TSNAXIP1 protein (Cat# [TP310670]). The protein was produced from HEK293T cells transfected with TSNAXIP1 cDNA clone (Cat# RC210670) using MegaTran 2.0 (Cat# [TT210002]).