

Product datasheet for **RC210684**

LOC149950 (NM_001010976) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC149950 (NM_001010976) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LOC149950
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC210684 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATTGTAATCCCCAAAGCGTGCAAACATCTGCCCCATCTCCCTGGAGGTGGCCTGTTCTCCACGC
CACCTCACAGGCTGCCTGGCGGACTCTCTCACAGCACTGTGCTCCCGGGGCCACTTGTACAGGCC
CATGCGAGAAGGACCTCGGGCAGTGTATAATCTCCAAGGGCCCACAGGAACAGCAGTGACAACCTGTGT
ATGAAGCACCTACTATGTGCTGGGACAAGAATGGAACAAGAAGACATGCTTTGCCCTCACCCCTGGAGG
GCAGCTTCCAGCCAGGGAGGCAAATACCTCTCCAGACTCTTCCACAGATCCCCAAACTCTGCCTCT
CTCCTTCAGATCCCTGCTCAGATGTCACCAGCTCTGTGCAGCCTCCCTGCCACCAAGTCTAAAATTGCCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC210684 protein sequence Red=Cloning site Green=Tags(s)
-------------------	--

MDCKSPKRANICPHLPGGGLFSTPPSQAAWRTLLTALCFPGPTCTGPMREGPRAYNPPRAHRNSSDNCV
MKHLLCAGDKNGTRRHLPSPLEGSFQPGRIPPPQTPSTDPQLPLSFRSLLRCHQLCAASLPPSLKLP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6544_g12.zip

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001010976

ORF Size: 420 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001010976.1](#), [NP_001010976.1](#)

RefSeq Size: 846 bp


RefSeq ORF: 422 bp

Locus ID: 149950


Cytogenetics: 20q11.21

MW: 15.1 kDa

Product images:



Circular map for RC210684



Western blot validation of overexpression lysate (Cat# [LY423249]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC210684 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).