

Product datasheet for **RC210701A1V**

Human C1orf187 (DRAXIN) (NM_198545) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human C1orf187 (DRAXIN) (NM_198545) AAV Particle
Tag:	Myc-DDK
Symbol:	C1orf187
Synonyms:	AGPA3119; C1orf187; neucrin; UNQ3119
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC210701 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTGGGCCTGCCATCCACACCGCTCCCATGCTGTTCTCTCGTCTCTGCTGCCCTGGAGCTGAGCC
TGGCAGGCGCCCTTGACCTGGGACCCCTGCCCGAACCTCCCTGAGAATCACATTGACCTCCCAGGCC
AGCGCTGTGGACGCCTCAGGCCAGCCACCACCGCGGGGGCCCGGCAAGAAGGAGTGGGGCCAGGC
CTGCCAGCCAGGCCAGGATGGGGCTGTGGTCACCGCCACAGGCAGGCCTCCAGGCTGCCAGAGGCTG
AGGGGCTGTGCTGAGCAGAGTCTGCAGGCCTGCTGCAGGACAAGGACCTGCTCTGGGACTGGCATT
GCCCTACCCCGAGAAGGAGAACCAGCTCCAGGTTGGGAGAGGACCAGGAAACGCAGCAGGGAGCACAAG
AGACGCAGGGACAGGTTGAGGCTGCACCAAGGCCGAGCCTTGGTCCGAGGTCCCAGCTCCCTGATGAAGA
AGGCAGAGCTCTCCGAAGCCCAGGTGCTGGATGCAGCCATGGAGGAATCCTCCACCAGCCTGGCGCCAC
CATGTTCTTTCTACCACCTTTGAGGCAGCACCTGCCACAGAAGAGTCCCTGATCCTGCCCGTCACCTCC
CTGCGGCCCCAGCAGGCACAGCCAGGTCTGACGGGGAGGTGATGCCACGCTGGACATGGCCTTGTTTCG
ACTGGACCGATTATGAAGACTTAAACCTGATGGTTGGCCCTCTGCAAAGAAGAAAGAGAAACACCGCGG
TAAACTCTCCAGTGATGGTAACGAAACATCACCAGCCGAAGGGGAACCATGCGACCATCACAAGACTGC
CTGCCAGGGACTTGCTGCGACCTGCGGGAGCATCTCTGCACACCCACAACCGAGGCCTCAACAACAAAT
GCTTCGATGACTGCATGTGTGTGAAGGGCTGCGCTGCTATGCCAAATTCACCGGAACCGCAGGGTTAC
ACGGAGGAAAGGCGCTGTGTGGAGCCGAGACGGCCAACGGCGACCGAGGATCCTTCATCAACGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210701 protein sequence Red=Cloning site Green=Tags(s) MAGPAIHTAPMLFLVLLLPLELSLAGALAPGTPARNLPENHIDLPGPALWTPQASHHRRRGPGKKEWGPGLPSQAQDGAVVTATRQASRLPEAEGLLPEQSPAGLLQDKDLLGLALPYPEKENRPPGWERTKRKRSREHKRRRDRLRLHQGRALVRGPSSLMKKAELSEAQVLDAAMEESSTSLAPTMFFLTTFEAPATEESLILPVTSLRPQQAQPRSDGEVMPITLDMALFDWTDYEDLKPDGWPSAKKKEKHRGKLSSDGNETSPAEGEPCDHHQDCLPGTCCDLREHLCTPHNRGLNNKCFDDCMCVEGLRCYAKFHRNRRVTRRKGRCEPETANGDQGSFINV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_198545
ORF Size:	1047 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_198545.3
RefSeq Size:	1773 bp
RefSeq ORF:	1050 bp
Locus ID:	374946
UniProt ID:	Q8NBI3
Cytogenetics:	1p36.22
MW:	38.6 kDa