

Product datasheet for **RC210721L4V**

WW domain binding protein 4 (WBP4) (NM_007187) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	WW domain binding protein 4 (WBP4) (NM_007187) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WW domain binding protein 4
Synonyms:	FBP21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007187
ORF Size:	1128 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210721).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007187.3
RefSeq Size:	2354 bp
RefSeq ORF:	1131 bp
Locus ID:	11193
UniProt ID:	O75554
Cytogenetics:	13q14.11
Domains:	WW, ZnF_U1


[View online »](#)

MW: 42.5 kDa

Gene Summary: This gene encodes WW domain-containing binding protein 4. The WW domain represents a small and compact globular structure that interacts with proline-rich ligands. This encoded protein is a general spliceosomal protein that may play a role in cross-intron bridging of U1 and U2 snRNPs in the spliceosomal complex A. [provided by RefSeq, Jul 2008]