

Product datasheet for **RC210745A1V**

Human SPCS2 (NM_014752) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SPCS2 (NM_014752) AAV Particle
Tag: Myc-DDK
Symbol: SPCS2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210745 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGGCGGCAGCTGTACAGGGCGGGAGAAGCGGTGGTAGCGGAGGCTGTAGTGGGCTGGTGGTGCTT
CCAACTGCGGGACAGGAAGTGGCGTAGCGCTTGTGGATAAGTGAAGATAGATGATAAGCCTGTAAA
AATTGACAAGTGGGATGGATCAGCTGTGAAAACTCTTTGGATGATTCTGCCAAAAGGTAAGTCTCTGAA
AAATACAAATATGTGGAGAATTTGGTCTAATTGATGGTCGCCTCACCATCTGTACAATCTCCTGTTCT
TTGCCATAGTGGCTTTGATTTGGGATTATATGCACCCCTTCCAGAGTCAAACCCGTTTGGCTTTGTG
TGTCATATCCTATTTTGTGATGATGGGGATTCTGACCATTATACCTCATATAAGGAGAAGAGCATCTTT
CTCGTGGCCACAGGAAGATCCTACAGGAATGGATCCTGATGATATTTGGCAGCTGTCCTCCAGTCTTA
AAAGGTTTGTGACAAATACACCTTGAAGCTGACCTTCATCAGTGGGAGAAAGCAGCAGCGGAAGC
CGAGTTCAAAAGTCCATTGCTAAGTTTTTGAACACAGTGGGACACTGGTCATGGATGCATATGAGCCT
GAAATATCCAGGCTCCATGACAGTCTTGCCATAGAAAGAAAAATAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC210745 protein sequence
Red=Cloning site Green=Tags(s)

MAAAVQGGRSGSGGCSGAGGASNCGTGSGRGLLDKWKIDDKPVKIDKWDGSAVKNSLDDSAKKVLL
KYKYVENFGLIDGRLTICTISCFFAIVALIWDYMHPPESKPVLALCVISYFVMMGILTIYTSYKEKSIF
LVAHRKDPTGMDPDDIWQLSSSLKRFDDKYTLKLTTFISGRKQQREAEFTKSIKFFDHSGTLVMDAYEP
EISRLHDSLAIERKIK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_014752
ORF Size:	678 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014752.2</u>
RefSeq Size:	2731 bp
RefSeq ORF:	681 bp
Locus ID:	9789
UniProt ID:	<u>Q15005</u>
Cytogenetics:	11q13.4
MW:	25 kDa